



24350 State Road 23 South, South Bend IN 46614

Tel: (574)233-4500 Fax: (574)288-4550

Mattresses: www.derbyindustries.com

Cases: www.derbycase.com

Email: sales @derbyindustries.com

Cleaning Procedures For The Derby Sealed Seam Mattress

We at Derby Industries **once again** strongly recommend the following cleaning procedures to ensure the durability of your quality Derby mattresses:

Use a hospital-grade cleaner or bleach-water solution on the entire mattress.

Follow with a final rinse of water only.

After cleaning, DRY THE MATTRESS COVER THOROUGHLY. Do not place a wet mattress into bed pans, as moisture may cause fabric damage over time.

NOTE: This damage IS NOT COVERED BY OUR WARRANTY.